



## Designing Antisense Oligonucleotides

### Contents


Introduction .....	1
Choosing a modified oligonucleotide .....	3
Ribonuclease H–Mediated Antisense Activity.....	3
Phosphorothioates .....	4
2' Substitutions .....	5
Target-Site Selection.....	7
Sequence Checks .....	9
Conclusions .....	10
References .....	11

### Introduction

Antisense oligonucleotides have been used for a number of years to modify the expression of specific genes both *in vivo* and *in vitro* [1]. The most potent mode of antisense activity is through RNase H–mediated degradation of RNA (Figure 1). The RNase H–endonuclease specifically cleaves RNA only when it is hybridized as a heteroduplex with DNA. In general, oligonucleotides that cannot invoke RNase H pathways are 10–100 fold less potent than RNase H–activating oligonucleotides of identical sequence.


Ideally, any researcher should be able to choose a specific target sequence of interest, order the synthesis of their designer antisense oligonucleotide, introduce it into their system of choice, and observe the effects. However, in practice, things are not this simple. The choice of sequence and chemistry of the antisense oligonucleotide is crucial to the success of the experiment [2].

**Problem:** Natural phosphodiester oligodeoxynucleotides (or oligoribonucleotides) are quickly digested by nucleases both *in vitro* and *in vivo*. Although multiple endo- and exonucleases exist and may be important *in vivo* [3], it appears that the bulk of biologically significant nucleolytic




activity in serum is seen as a 3' exonuclease [4] (Figure 2). Within the cell nucleolytic activity is seen as both 5' and 5' exonuclease [5] (Figure 3).

Degradation of Antisense Oligonucleotides in Plasma Occurs Solely Through a 3' - exonuclease Activity



**Figure 2. Substrate specificity and kinetics of degradation of antisense oligonucleotides by 3' exonuclease in plasma.** 1=control, 2=15 minutes, 3=30 minutes, 4=1 hour, 5=2 hours, \*=32P, ^=TCDM phosphotriester. Adapted from [4].

Degradation of Antisense Oligonucleotides in Cells Occurs by Both 3' and 5' Exonuclease Activities



**Figure 3. Pathways of degradation and mechanism of action of antisense oligonucleotides in *Xenopus laevis* embryos.** 1 = control, 2 = 1 minute, 3 = 10 minutes, 4 = 20 minutes, 5 = 40 minutes, 6 = 60 minutes. Oligonucleotides modified at both the 3' and 5' ends by methoxyethylphosphoramidate are protected from exonuclease digestion. Oligonucleotides modified at only the 3' end are digested by a 5' exonuclease. Adapted from [5].

**Solution:** Synthetic antisense oligonucleotide probes differ from natural DNA or RNA in that they are modified to some form that is nuclease-resistant. Many such modifications have been conceived and tested, few of which work as well as intended. Ideally, these compounds should:

- Activate RNase H–degradation pathways.
- Be easy to make and inexpensive.
- Not be physiologically toxic.
- Not be easily degraded.
- Not disrupt normal Watson–Crick base-pairing.
- Not induce any unanticipated sequence-independent biologic effects.

## Choosing a modified oligonucleotide

Many different modifications have been substituted by investigators into the native phosphodiester oligodeoxynucleotide polymer to limit nuclease sensitivity in an effort to enhance their utility as antisense agents. Short oligos have very brief intracellular half-lives regardless of modifications, but modified oligos with a length of 16 or so residues can have half-lives on the order of days [3]. Residue modifications can also affect the ability of an oligo to trigger RNase H mediated degradation of RNA following hybrid formation. As RNase H mediated destruction of RNA may be the primary mode of antisense action, this effect cannot be ignored. Human RNase H cleaves RNA in RNA:DNA hybrids. The enzyme requires divalent cations to be active. Peak activity occurs in the presence of 10mM Mg<sup>++</sup>, 5mM Co<sup>++</sup>, or 0.5mM Mn<sup>++</sup>. Action of the enzyme leaves a 5'-phosphate and a 3'-hydroxyl [4].

### Summary of Candidate Compounds for Potential Use as Antisense Oligonucleotide Probes

---

Chemistry	Resistance To Nucleases	Activation of RNase H pathways	Comments
Phosphodiester	NO	YES	Natural DNA
Methylphosphonate.	YES	NO	Will work as chimera
Phosphorothioate	YES	YES	Stimulates many biologic responses in a sequence non-specific fashion
$\alpha$ -nucleoside	YES	NO	Will work as $\alpha\beta\alpha$ chimera
2'-O-substituted RNA	YES	NO	Greater stability than DNA:RNA hybrids Will work as chimera
Phosphoramidite	YES	NO	Will work as chimera
Morpholino	YES	NO	Will work as steric blocker
Chimeras	YES	YES	Contain an internal core of unmodified phosphodiester DNA flanked by modified residues

## Ribonuclease H–Mediated Antisense Activity


The RNase H class of endonucleases acts primarily in the nucleus, although activity can be detected in cytoplasm [6]. It is thought that the antisense oligonucleotide probe binds specifically to target mRNA, which initiates RNase H–mediated degradation of the double-

stranded antisense probe:mRNA hybrid [5] (Figure 5). RNA degradation products corresponding to the fragments expected from RNase H action can be detected in living cells treated with antisense agents [7]. The precise mechanism of action is not known but it has been proposed to be similar to that of DNase I [8]. Most DNAs that have been modified to be nuclease-resistant do not form a heteroduplex structure. One exception is the phosphorothioate-modified oligonucleotides discussed below.

Human RNase H enzymes will completely digest RNA in heteroduplex form and have optimal activity when the gap between phosphorothioate linked bases and phosphodiester linked bases is at least 6–8 residues [9]. However, when enzyme levels are limiting, it becomes clear that certain sites cleave preferentially [9, 10, 11, 12]. Cleavage rates of a single ribonucleotide residue embedded in DNA proceeds at A > U > C > G; the rate of cleavage at rA is over four times faster than at an rG [4] (Figure 3). Further, the RNA secondary and tertiary structures may present certain sites as topologically favored for enzymatic attack. The scientist must be certain that the target sequence is accessible within the folded RNA structure when designing antisense compounds.

## Phosphorothioates

The easy-to-synthesize phosphorothioate oligonucleotides (Figure 4) assume the native Watson–Crick nucleotide hydrogen-bonding patterns, can activate RNase H-mediated degradation of cellular mRNA, and are nuclease-resistant [13, 14, 15]. The antisense effects of the phosphorothioates can be observed for over 48 hours after a single application to tissue culture cells [16]. This degree of stability is needed for *in vivo* work [17]. However, the actual stability of a phosphoro-thioate oligonucleotide in a specific experiment can vary with each sequence and cell line examined [18]. Although early work using these compounds was very encouraging, it has become clear that some of the most exciting results were actually due to sequence independent biological effects of phosphorothioate DNA (sulfated polyanion) and did not result from true antisense mechanisms [19, 20].



**Figure 4. Phosphorothioate Oligonucleotide (GC)**

The synthesis of chimeric oligos having modified 3'- and 5'- ends with a normal phosphodiester central core has been employed using phosphoramidites and methylphosphonates with good results. Similar chimeras have been evaluated for phosphorothioate-modified oligos in an attempt to lower the total thiol-modified content of the oligo and decrease side-effects/toxicity. Leaving a central core of as few as 3 phosphodiester linkages leaves the oligo relatively unprotected from *Bal31* nuclease digestion [21]. Such end-blocked oligos may still be protected sufficiently to be useful in practice, however, as the predominant activity that destroys oligos in serum is a 3'-exonuclease. Chimeric structure has been shown to have decreased side effect profile *in vivo* [22, 23]. In spite of some disadvantages, phosphorothioates remain the most popular agent used in antisense studies. When they work, they work well and are cost effective


The phosphorothioate modification slightly decreases the relative binding affinity of an oligo, similar to most other phosphate modifiers, such as phosphoramidites or methylphosphonates. Other changes can be introduced into such oligos to increase their relative binding affinities and compensate for this decrease in  $T_m$ . The addition of C-5 propyne pyrimidines into phosphorothioate oligos can increase their binding affinity and thereby their potency such that they become effective in the 100nM range. This shift in the dose response curves might lower the effective dosage of these agents to the point that worries over toxic, non-specific side effects will disappear [16, 24, 25]. The use of C-5 propyne-modified oligos allow for a decrease in oligo length. Oligos as short as 11 bases can have potent antisense effects [26]. A similar enhancement can be achieved using C-7 modified purines - 7-deaza-2'-deoxyguanosine and 7-deaza-2'-deoxyadenosine [27]. Further improvement in oligo  $T_m$  and nuclease stability can be achieved through the combined use of both 2'-*O*-modification plus thiolation [28].

Like all mono-substituted phosphate backbone modifiers, phosphorothioates are chiral. Chirally pure phosphorothioates can be synthesized. Pure *R<sub>p</sub>* configuration oligos have higher binding affinity (higher  $T_m$ ) than pure *S<sub>p</sub>* or mixed congeners. The *R<sub>p</sub>* species was found to have increased RNase-H-activating potential than *S<sub>p</sub>* or chirally mixed oligos [29].

## 2' Substitutions

Substitutions at the ribose 2' position can significantly alter an oligo's nuclease stability and binding properties [30]. Hundreds of such modifications have been synthesized and tested for activity, including 2'-fluoro, 2'-*O*-methyl, 2'-*O*-allyl, 2'-*O*-propyl, and 2'-*O*-pentyl. These agents confer varying degrees of nuclease stability to an oligo. 2'-*O*-allyl analogs have a life span in serum of about 24 hours. 2'-*O*-methyl analogs (Figure 5) are some 10 fold less stable than the allyl class but are nonetheless more stable than native phosphodiester oligos [3, 31]. This degree of stability is not sufficient for most antisense applications. Oligos of this class do not activate RNase H [32], but can have antisense effects through a steric blocking mode of action. 2'-*O*-alkyl


oligoribonucleotides have unambiguously been shown to be capable of inhibiting mRNA translation through steric effects *in vitro* [33]. The introduction of substitutions that increase binding affinity can enhance steric antisense effects. A C-5 propyne-pyrimidine 2'-*O*-allyl oligo was described which had a  $T_m$  of  $>90^\circ\text{C}$  for its target RNA. This construct showed antisense effect after microinjection into CV1 cells with an  $\text{IC}_{50}$  of 5  $\mu\text{M}$ . However, an RNase H-activating C-5 propyne phosphorothioate oligo of identical sequence had a  $T_m$  of only  $79^\circ\text{C}$  but showed an  $\text{IC}_{50}$  of 0.25  $\mu\text{M}$ , a 20 fold increase in antisense potency [25]. RNase H activating capacity can be introduced into 2'-substituted oligos through the construction of 2'-modified / 2'-deoxy gap chimeras [34]. These chimeras can also be thioated to enhance nuclease resistance [28].



**Figure 5. 2'-*O*-methyl oligonucleotide (AG)**

Oligo-2'-fluoro-2'-deoxynucleotide phosphoramidites have been synthesized and these oligos appear to offer the best increase of  $T_m$  for this class of compounds [35]. Unfortunately, the addition of bulkier groups at the 2'-position is needed to confer nuclease resistance but these bulkier groups provide less benefit in improved binding affinities. The relative binding affinity of these compounds is: fluoro > methoxy > propoxy > pentoxy = deoxy. Conversely, the relative nuclease resistance of these compounds is: pentoxy > propoxy > methoxy > fluoro = deoxy [21]. Engineering a self-complementary hairpin onto the 3' -end can increase nuclease resistance of these compounds [36].

Griffey et al. [37] report the synthesis of a 2'-O-aminopropyl modified oligos. This substitution was more nuclease resistant than phosphorothioates and had slightly increased binding affinity compared to the standard phosphodiester oligo. "Gapmers", incorporating 2'-modified flanks with a phosphorothioate core, were synthesized and tested for antisense activity. The new constructs had 5-10 fold more activity (in tissue culture) than an otherwise identical homogenous phosphorothioate compound. One unusual 2'-substitution has recently been described. (2'S)-2'-deoxy-2'-C-methyl oligonucleotides are resistant to *Bal31*, P1, and mung bean nucleases and are stable in serum [38]. Relative binding affinities and potential for RNase H activation are not known. These compounds are now considered the reagents of choice for "second generation" compounds (Figure 6).



**Figure 6. Recommended Chimeric Antisense Oligonucleotide Design**


Phosphorothioates show increased binding to cellular proteins and components of the extracellular matrix as compared to natural phosphodiester oligonucleotides [19, 20]. This binding appears to be due to the polyanionic nature of these compounds; they behave similar to dextran sulfate and heparin sulfate. This binding can displace or mimic the binding of natural ligands to assorted proteins, such as receptors or adhesion molecules [39]. In fact, any of the heparin-binding class of proteins may also bind phosphorothioates [40, 41]. Phosphodiester DNA is a polyanion and may nonspecifically bind proteins, but due to nuclease action has such a shortened lifespan that the impact of this effect is most likely limited.

Chimeric oligonucleotides have been synthesized with phosphoramidite-modified 3'- and 5'-ends with a normal phosphodiester central core (Figure 7). Leaving a central core of as few as 3 phosphodiester linkages leaves the oligonucleotide relatively unprotected from *Bal31* nuclease digestion [21]. However, since the predominant activity that destroys

oligonucleotides in serum is a 3'-exonuclease, end-blocked oligonucleotides are protected sufficiently to be useful in practice. Chimeric structure has been shown to decrease side effects *in vivo* [22, 23].

### Evolution of Modern Chimera Design Gap Length

1 2 3 4 5 6 7 8



Lane	Oligo		Gap	Degradation
1-4	Controls		-	No
5	An2-8	C*T*G*A <sub>0</sub> C <sub>0</sub> A <sub>0</sub> A <sub>0</sub> C <sub>0</sub> A <sub>0</sub> T <sub>0</sub> G <sub>0</sub> A*C*T*G*C	3-8-4	Yes
6	An2-6	C*T*G*A <sub>0</sub> C <sub>0</sub> A <sub>0</sub> A <sub>0</sub> C <sub>0</sub> A <sub>0</sub> T <sub>0</sub> G*A*C*T*G*C	4-6-5	Yes
7	An2-5	C*T*G*A*C <sub>0</sub> A <sub>0</sub> A <sub>0</sub> C <sub>0</sub> A <sub>0</sub> T <sub>0</sub> G*A*C*T*G*C	5-5-5	No
8	An2-4	C*T*G*A*C*A <sub>0</sub> A <sub>0</sub> C <sub>0</sub> A <sub>0</sub> T <sub>0</sub> G*A*C*T*G*C	5-4-6	No

**Figure 7. Pathways of degradation and mechanism of action of antisense oligonucleotides in *Xenopus Laevis* embryos.** \* Phosphorothioate internucleotide linkage, 0 phosphodiester internucleotide linkage. Central core of 6-8 phosphodiester linkages allow degradation by RNase H (Adapted from [5]).

### Target-Site Selection

Finding the right target site may be easy or the most difficult step in the whole antisense process. It ultimately determines the success of any antisense investigation. This section discusses some of the things that can cause problems. Lest these appear daunting or even overwhelming, it should be pointed out that scientists at companies such as Isis and Sequitur claim that they achieve close to 100% success in finding useful antisense oligonucleotides for given targets.

Some investigators have found that about one out of ten phosphorothioate 20-mers will exhibit good activity and that this can be improved using superior chemical modifications [42]. If the oligomers are randomly selected, only a small fraction may be active [43, 44]. Sometimes only a single oligonucleotide, most often one complementary to the initiator AUG, will inhibit expression of a particular RNA. Obviously, testing as many oligonucleotides as possible improves the chance of identifying one that is particularly potent. This can be expedited by monitoring gene expression using RT-PCR in a 96-well format or some other high-throughput procedure such as micro-array technology [42].

**Length.** The stability of hybrids can depend on length, particularly for shorter oligonucleotides. Also, as the length of an oligonucleotide increases, it is less likely that it will encounter a complementary sequence other than the targeted RNA. On the other hand, increasing the length of the oligonucleotide increases the probability that it will bind a partially complementary sequence in a nontarget message, thereby activating RNase H, which requires only six or seven base pairs in a heteroduplex substrate for activation. The usual length for antisense oligonucleotides is around 20 bases, which is a convenient size for synthesis and long enough, on statistical grounds, to be unique in the human genome. In special cases, structural features in the target RNA may enable the use of shorter oligomers [24, 45]. However, most studies find a decrease or loss of antisense activity as length is reduced from twenty to ten bases [46].

**Conformational and thermodynamic considerations.** The major problem lies with the secondary and tertiary folding that can make much of the RNA inaccessible to a molecule as large as an oligonucleotide. Even those sequences that appear to be accessible may already be involved in intramolecular hydrogen bonding, stacking interactions, or in solvation that would be disrupted by hybridization of an oligonucleotide. Consequently, hybridization-induced rearrangement of the existing RNA structure may carry a prohibitive thermodynamic penalty. On the other hand, single-stranded sequences within the RNA may be preordered by stacking into helical conformations that are particularly favorable for hybridization. The exceptional stability of hybrids formed between the loops of two hairpins (kissing interactions) is well known and is important in the association of natural antisense RNAs with their targets. Factors such as these may account for the preferential hybridization of oligonucleotides to the 5' rather than 3' side of loops [47, 48, 49]. The difficulty of predicting accessible binding sites is illustrated by a study comparing experimentally determined RNase H-accessible sites in four RNA species with those predicted by known secondary or tertiary structure [50]. Overall, most cleavage reactions occurred at single-stranded sites in the RNA. While some single-stranded sites were good targets, many were not, and many good sites were located in double-stranded regions.

Even though the rules for base-pairing are very simple, additional subtleties govern the hybridization of oligonucleotides to RNA that are not well understood. The behavior of oligonucleotides is very dependent on the terminal nucleotides. Moreover, small changes in the length or a shift in binding site of one or two nucleotides can profoundly affect the kinetics of hybrid formation [51, 52, 53, 54, 55]. Even a few base changes that do not change the thermodynamic stability of the duplex may greatly change the kinetics of hybridization [56, 57, 58]. These effects may partially account for the efficacy of different antisense oligonucleotides *in vivo*.

Sophisticated modeling programs such as Foldsplit [59, 60] can often yield more active antisense oligonucleotides than those obtained by random selection [58]. Such modeling programs have predicted that, on average, there is about one good target site per 1000 bases in RNA.

**Targeting protein-binding sites.** Sequences in RNA that interact with proteins, ribosomes, spliceosomes, and other large entities are also likely to be accessible to oligonucleotides,

assuming no unwinding activity is required. Early on, the cap, initiator codon, and 3'-end were selected as targets [61]. Many later studies have also found that the initiator codon is a good target and has become something of an industry standard despite the occasional failure. Antisense oligonucleotides have also been used to redirect splicing to prevent formation of the functional mRNA [62]. Successful targeting of splice sites requires that the oligonucleotide gain access to the nucleus whereas inhibition of translation may be accomplished by hybridization in the cytoplasm. In some cases, it is possible to use oligomers that do not induce cleavage, particularly if they hybridize strongly such as the morpholino analogues or 2'-O-methyl derivatives [63, 64, 65, 66].


**Testing libraries of oligomers.** cDNA libraries may be made more efficacious for screening by cutting the cDNA into fragments of approximately ten bases [50]. RNase H is often used to identify those sequences showing greatest activity [50, 56, 57, 67, 68]. These studies can yield both kinetic and equilibrium hybridization data and have been shown to be predictive of activity in cells [56, 68, 69]. These experiments can also be performed in cell extracts to more closely mimic the structure and protein binding of the RNA in cells [70]. A potential problem using libraries of random oligonucleotides is that members of the library may interact with each other inhibiting hybridization with the target [47]. Alternatively, they could assist each other to hybridize to the target through cooperative binding [71, 72, 73]. This problem can be avoided by using arrays of immobilized oligonucleotides [53].

### Sequence Checks

Before synthesizing an antisense oligonucleotide, the investigator should check the sequence for various features that could affect its activity. For instance, if the sequence complements nontarget RNAs, the probe may not be useful. In addition, the oligonucleotide should be examined for self-complementary sequences that might interfere with hybridization to the target. Certain sequence motifs have potent biological effects unrelated to antisense activity. This was unrecognized by early investigators and is still the biggest pitfall for the unwary researcher.

**CG-containing sequences.** Unmethylated CG sequences occur more frequently in bacterial DNA than they do in eukaryotic DNA. Consequently, their presence may be used as a signal of bacterial infection by the immune system [74]. Oligonucleotides containing CG can act as immunostimulators by causing proliferation of B lymphocytes; by activating macrophages, dendritic cells, and T cells; and by inducing cytokine release [20, 74, 75, 76, 77, 78]. These CG-mediated immune effects depend on the sequences flanking the CG dimer, and are strongest with the purine.purine.CG.pyrimidine.pyrimidine motif [20]. These CG effects occur with phosphorothioates as well as with phosphodiester [79], and may be responsible for some of the activities of oligonucleotides reported *in vivo*. The easiest solution is to choose oligonucleotides that do not contain CG, particularly those with flanking sequences that favor immune stimulation. An alternative is to replace the C in CG sequences with 5-methylcytidine [80]. Although it increases the expense of synthesis, this 5-methylcytidine substitution prevents immune stimulation without affecting hybridization.

**Tetraplexes.** In addition to forming duplexes and hairpins by Watson–Crick base pairing, some oligonucleotides can form structures comprising three, four, or more strands. In particular, formation of tetraplexes with potent biological activity has caused some problems in the antisense field [81, 82]. Investigators should carefully examine all oligonucleotides very rich in a particular nucleoside, particularly if they show repeated sequences or have multiple occurrences of two or more adjacent identical bases. Oligomers with multiple repeats of two or more consecutive Gs or Cs may form tetraplexes and other non–Watson–Crick structures. Not all oligomers with such features will necessarily form these higher order structures, particularly in physiological conditions. Nonetheless, such sequences raise warning flags and there is a well-documented danger in ascribing biological effects to an antisense mechanism without careful investigation.



**Figure 8. Possible Guanine Tetraplex Structure**

The most extensively studied tetraplexes are formed by oligonucleotides containing multiple adjacent guanine residues (Figure 8). These may occur in a single run of around four residues but they can also be found in repeated GG or GGG motifs that occur in close proximity [83, 84]. Even if they do not form tetraplexes, G-rich sequences with multiple GG dimers may form other unusual structures depending on sequence context [85]. Tetraplex-forming runs of Gs seem to have an affinity for various proteins and when included in synthetic oligonucleotides, they produce a multitude of biological effects. For example, researchers have identified tetraplexes that bind to thrombin [86] and to the HIV envelope protein [87]. Other tetraplexes have been shown to bind to transcription factors [88] or to produce antiproliferative effects by protein binding [89]. The ability to form tetraplexes can be blocked by replacing guanosine residues with 7-deazaG [90] or 6-thioG [91]. It should also be noted that a phosphorothioate oligonucleotide containing only C residues was shown to have activity similar to one containing a G-tetraplex [82].

**Other motifs.** Investigators have suggested that stretches of purines in the target might stabilize the heteroduplex formed [92]. From examining the the sequences of active antisense oligonucleotides in many published studies, investigators have proposed that selecting a target containing the sequence GGGA gives a much better chance of success [93].

## Conclusions

Antisense oligonucleotides can provide a tremendous amount of fundamental data to committed investigators. Systematic down regulation of particular gene products can be key in determining a gene's biological function and its role in a metabolic or regulatory pathway. Similarly, antisense technology provides a valuable tool for discerning genetic contributions to disease. The potential for antisense-based oligonucleotide drugs is only now beginning to be realized. Isis Pharmaceuticals recently brought to market Vitravene™, an antisense oligonucleotide drug for CMV retinitis. A number of other antisense oligonucleotides are currently in clinical trials, including those for treating malignancies and for targeting diseases such as psoriasis, hepatitis C, and rheumatoid arthritis. It is only through the proper design of

antisense oligonucleotides that these useful tools have come to yield a plethora of meaningful results, both in pure and applied research.

## References

1. Scanlon KJ, Ohta Y, et al. (1995) Oligonucleotide-mediated modulation of mammalian gene expression. *FASEB J*, 9: 1288–1296.
2. Stein CA and Cheng YC. (1993) Antisense oligonucleotides as therapeutic agents — is the bullet really magical? *Science*, 261: 1004–1012.
3. Fisher TL, Terhorst T, et al. (1993) Intracellular disposition and metabolism of fluorescently-labeled unmodified and modified oligonucleotides microinjected into mammalian cells. *Nucleic Acids Res*, 21: 3857–3865.
4. Eder PS, DeVine RJ, et al. (1991) Substrate specificity and kinetics of degradation of antisense oligonucleotides by a 3' exonuclease in plasma. *Antisense Res Dev*, 1: 141–151.
5. Dagle JM, Weeks DL, and Walder JA. (1991) Pathways of degradation and mechanism of action of antisense oligonucleotides in *Xenopus laevis* embryos. *Antisense Res Dev*, 1: 11–20.
6. Cazenave C, Frank P, et al. (1994) Characterization and subcellular localization of ribonuclease H activities from *Xenopus laevis* oocytes. *J Biol Chem*, 269: 25185–25192.
7. Giles RV, Spiller DG, and Tidd DM. (1995) Detection of ribonuclease H-generated mRNA fragments in human leukemia cells following reversible membrane permeabilization in the presence of antisense oligodeoxynucleotides. *Antisense Res Dev*, 5: 23–31.
8. Ishikawa K, Okumura M, et al. (1993) Crystal structure of ribonuclease H from *Thermus thermophilus* HB8 refined at 2.8 Å resolution. *J Mol Biol*, 230: 529–542.
9. Eder PS and Walder JA. (1991) Ribonuclease H from K562 human erythroleukemia cells. Purification, characterization, and substrate specificity. *J Biol Chem*, 266: 6472–6479.
10. Eder PS, Walder RY, and Walder JA. (1993) Substrate specificity of human RNase H1 and its role in excision repair of ribose residues misincorporated in DNA. *Biochimie*, 75: 123–126
11. Toulme JJ, Boiziau C, et al. (1996) Ribonuclease H: from enzymes to antisense effects of oligonucleotides. In *DNA and RNA Cleavers and Chemotherapy of Cancer and Viral Diseases*, Meunier B, ed. Amsterdam, Netherlands: Kluwer Academic Publishers, 271–288.
12. Frank P, Albert S, et al. (1994) Purification and characterization of human ribonuclease HIII. *Nucleic Acids Res*, 22: 5247–5254.
13. Stein CA, Tonkinson JL, and Yakubov L. (1991) Phosphorothioate oligodeoxynucleotides — anti-sense inhibitors of gene expression? *Pharmacol Ther*, 52: 365–384.
14. Crooke ST and Lebleu B, eds. (1993) *Antisense Research and Applications*. Boca Raton FL: CRC Press.
15. Srinivasan SK and Iversen P. (1995) Review of in vivo pharmacokinetics and toxicology of phosphorothioate oligonucleotides. *J Clin Lab Anal*, 9: 129–137.

16. Bonham MA, Brown S, et al. (1995) An assessment of the antisense properties of RNase H-competent and steric-blocking oligomers. *Nucleic Acids Res*, 23: 1197–1203.
17. Yaida Y and Nowak TS. (1995) Distribution of phosphodiester and phosphorothioate oligonucleotides in rat brain after intraventricular and intrahippocampal administration determined by in situ hybridization. *Regul Pept*, 59: 193–199.
18. Crooke ST, Lemonidis KM, et al. (1995) Kinetic characteristics of Escherichia coli RNase H1: cleavage of various antisense oligonucleotide-RNA duplexes. *Biochem J*, 312 ( Pt 2): 599–608.
19. Gura T. (1995) Antisense has growing pains. *Science*, 270: 575–577.
20. Krieg AM and Stein CA. (1995) Phosphorothioate oligodeoxynucleotides: antisense or anti-protein? *Antisense Res Dev*, 5: 241.
21. Monia BP, Johnston JF, et al. (1996) Nuclease resistance and antisense activity of modified oligonucleotides targeted to Ha-ras. *J Biol Chem*, 271: 14533–14540.
22. Retaux S, McNeill L, and Harris WA. (1996) Engrailed, retinotectal targeting, and axonal patterning in the midbrain during Xenopus development: an antisense study. *Neuron*, 16: 63–75.
23. Hebb MO and Robertson HA. (1997) End-capped antisense oligodeoxynucleotides effectively inhibit gene expression in vivo and offer a low-toxicity alternative to fully modified phosphorothioate oligodeoxynucleotides. *Brain Res Mol Brain Res*, 47: 223–228.
24. Wagner RW, Matteucci MD, et al. (1996) Potent and selective inhibition of gene expression by an antisense heptanucleotide. *Nat Biotechnol*, 14: 840–844.
25. Moulds C, Lewis JG, et al. (1995) Site and mechanism of antisense inhibition by C-5 propyne oligonucleotides. *Biochemistry*, 34: 5044–5053.
26. Flanagan WM, Kothavale A, and Wagner RW. (1996) Effects of oligonucleotide length, mismatches and mRNA levels on C-5 propyne-modified antisense potency. *Nucleic Acids Res*, 24: 2936–2941.
27. Buhr CA, Wagner RW, et al. (1996) Oligodeoxynucleotides containing C-7 propyne analogs of 7-deaza-2'-deoxyguanosine and 7-deaza-2'-deoxyadenosine. *Nucleic Acids Res*, 24: 2974–2980.
28. McKay RA, Cummins LL, et al. (1996) Enhanced activity of an antisense oligonucleotide targeting murine protein kinase C-alpha by the incorporation of 2'-O-propyl modifications. *Nucleic Acids Res*, 24: 411–417.
29. Koziolkiewicz M, Krakowiak A, et al. (1995) Stereodifferentiation--the effect of P chirality of oligo(nucleoside phosphorothioates) on the activity of bacterial RNase H. *Nucleic Acids Res*, 23: 5000–5005.
30. Lamm GM, Biencowe BJ, et al. (1991) Antisense probes containing 2-aminoadenosine allow efficient depletion of U5 snRNP from HeLa splicing extracts. *Nucleic Acids Res*. 19: 3193–3198.
31. Cummins LL, Owens SR, et al. (1995) Characterization of fully 2'-modified oligoribonucleotide hetero- and homoduplex hybridization and nuclease sensitivity. *Nucleic Acids Res*, 23: 2019–2024.

32. Lamond AI and Sproat BS. (1993) Antisense oligonucleotides made of 2'-O-alkylRNA: their properties and applications in RNA biochemistry. *FEBS Lett*, 325: 123–127.
33. Johansson HE, Belsham GJ, et al. (1994) Target-specific arrest of mRNA translation by antisense 2'-O-alkyloligonucleotides. *Nucleic Acids Res*, 22: 4591–4598.
34. Monia BP, Lesnik EA, et al. (1993) Evaluation of 2'-modified oligonucleotides containing 2'-deoxy gaps as antisense inhibitors of gene expression. *J Biol Chem*, 268: 14514–14522.
35. Schultz RG and Gryaznov SM. (1996) Oligo-2'-fluoro-2'-deoxynucleotide N3'-->P5' phosphoramidates: synthesis and properties. *Nucleic Acids Res*, 24: 2966–2973.
36. Kuwasaki T, Hosono K, et al. (1996) Hairpin antisense oligonucleotides containing 2'-methoxynucleosides with base-pairing in the stem region at the 3'-end: penetration, localization, and Anti-HIV activity. *Biochem Biophys Res Commun*, 228: 623–631.
37. Griffey RH, Monia BP, et al. (1996) 2'-O-aminopropyl ribonucleotides: a zwitterionic modification that enhances the exonuclease resistance and biological activity of antisense oligonucleotides. *J Med Chem*, 39: 5100–5109.
38. Iribarren AM, Cicero DO, and Neuner PJ. (1994) Resistance to degradation by nucleases of (2'S)-2'-deoxy-2'-C-methyloligonucleotides, novel potential antisense probes. *Antisense Res Dev*, 4: 95–98.
39. Khaled Z, Benimetskaya L, et al. (1996) Multiple mechanisms may contribute to the cellular anti-adhesive effects of phosphorothioate oligodeoxynucleotides. *Nucleic Acids Res*, 24: 737–745.
40. Guvakova MA, Yakubov LA, et al. (1995) Phosphorothioate oligodeoxynucleotides bind to basic fibroblast growth factor, inhibit its binding to cell surface receptors, and remove it from low affinity binding sites on extracellular matrix. *J Biol Chem*, 270: 2620–2627.
41. Stein CA. (1996) Phosphorothioate antisense oligodeoxynucleotides: questions of specificity. *Trends Biotechnol*, 14: 147–149.
42. Bennett CF and Cowser LM. (1999) Antisense oligonucleotides as a tool for gene functionalization and target validation. *Biochim Biophys Acta*, 1489: 19–30.
43. Peyman A, Helsenberg M, et al. (1995) Inhibition of viral growth by antisense oligonucleotides directed against the IE110 and the UL30 mRNA of herpes simplex virus type-1. *Biol Chem Hoppe Seyler*, 376: 195–198.
44. Monia BP, Sasmor H, et al. (1996) Sequence-specific antitumor activity of a phosphorothioate oligodeoxyribonucleotide targeted to human C-raf kinase supports an antisense mechanism of action in vivo. *Proc Natl Acad Sci U S A*, 93: 15481–15484.
45. Disney MD, Testa SM, and Turner DH. (2000) Targeting a *Pneumocystis carinii* group I intron with methylphosphonate oligonucleotides: backbone charge is not required for binding or reactivity. *Biochemistry*, 39: 6991–7000.
46. Cowser LM, Fox MC, et al. (1993) In vitro evaluation of phosphorothioate oligonucleotides targeted to the E2 mRNA of papillomavirus: potential treatment for genital warts. *Antimicrob Agents Chemother*, 37: 171–177.

47. Bruice TW and Lima WF. (1997) Control of complexity constraints on combinatorial screening for preferred oligonucleotide hybridization sites on structured RNA. *Biochemistry*, 36: 5004–5019.
48. Lima WF, Monia BP, et al. (1992) Implication of RNA structure on antisense oligonucleotide hybridization kinetics. *Biochemistry*, 31: 12055–12061.
49. Puri N and Chattopadhyaya J. (1999) How kinetically accessible is an RNA target for hybridization with an antisense oligo? A lesson from an RNA target which is as small as a 20mer. *Tetrahedron*, 55: 1505–1516.
50. Matveeva O, Felden B, et al. (1997) A rapid in vitro method for obtaining RNA accessibility patterns for complementary DNA probes: correlation with an intracellular pattern and known RNA structures. *Nucleic Acids Res*, 25: 5010–5016.
51. Rittner K, Burmester C, and Sczakier G. (1993) In vitro selection of fast-hybridizing and effective antisense RNAs directed against the human immunodeficiency virus type 1. *Nucleic Acids Res*, 21: 1381–1387.
52. Kronenwett R, Haas R, and Sczakier G. (1996) Kinetic selectivity of complementary nucleic acids: bcr-abl-directed antisense RNA and ribozymes. *J Mol Biol*, 259: 632–644.
53. Milner N, Mir KU, and Southern EM. (1997) Selecting effective antisense reagents on combinatorial oligonucleotide arrays. *Nat. Biotechnol*, 15: 537–541.
54. Venturini F, Braspenning J, et al. (1999) Kinetic selection of HPV 16 E6/E7-directed antisense nucleic acids: anti-proliferative effects on HPV 16-transformed cells. *Nucleic Acids Res*, 27: 1585–1592.
55. Sczakiel G. (2000) Theoretical and experimental approaches to design effective antisense oligonucleotides. *Front Biosci*, 5: D194–D201.
56. Lima WF, Brown-Driver V, et al. (1997) Combinatorial screening and rational optimization for hybridization to folded hepatitis C virus RNA of oligonucleotides with biological antisense activity. *J Biol Chem*, 272: 626–638.
57. Ho SP, Bao Y, et al. (1998) Mapping of RNA accessible sites for antisense experiments with oligonucleotide libraries. *Nat Biotechnol*, 16: 59–63.
58. Patzel V, Steidl U, et al. (1999) A theoretical approach to select effective antisense oligodeoxyribonucleotides at high statistical probability. *Nucleic Acids Res*, 27: 4328–4334.
59. Sczakiel G, Homann M, and Rittner K. (1993) Computer-aided search for effective antisense RNA target sequences of the human immunodeficiency virus type 1. *Antisense Res Dev*, 3: 45–52.
60. Sczakiel G and Tabler M. (1997) Computer-aided calculation of the local folding potential of target RNA and its use for ribozyme design. *Methods Mol Biol*, 74: 11–15.
61. Goodchild J, Carroll E 3<sup>rd</sup>, and Greenberg JR. (1988) Inhibition of rabbit beta-globin synthesis by complementary oligonucleotides: identification of mRNA sites sensitive to inhibition. *Arch Biochem Biophys*, 263: 401–409.

62. Kole R and Sazani P. (2001) Antisense effects in the cell nucleus: modification of splicing. *Curr Opin Mol Ther*, 3: 229–234.
63. Arora V, Knapp DC, et al. (2000) c-Myc antisense limits rat liver regeneration and indicates role for c-Myc in regulating cytochrome P-450 3A activity. *J Pharmacol Exp Ther*, 292: 921–928.
64. Iversen PL. (2001) Phosphorodiamidate morpholino oligomers: favorable properties for sequence-specific gene inactivation. *Curr Opin Mol Ther*, 3: 235–238.
65. Lacerra G, Sierakowska H, et al. (2000) Restoration of hemoglobin A synthesis in erythroid cells from peripheral blood of thalassemic patients. *Proc Natl Acad Sci U S A*, 97: 9591–9596.
66. Mercatante DR, Bortner CD, et al. (2001) Modification of alternative splicing of Bcl-x pre-mRNA in prostate and breast cancer cells. analysis of apoptosis and cell death. *J Biol Chem*, 276: 16411–16417.
67. Southern EM, Case-Green SC, et al. (1994) Arrays of complementary oligonucleotides for analysing the hybridisation behaviour of nucleic acids. *Nucleic Acids Res*, 22: 1368–1373.
68. Birikh KR, Berlin YA, et al. (1997) Probing accessible sites for ribozymes on human acetylcholinesterase RNA. *RNA*, 3: 429–437.
69. Ho SP, Britton DH, et al. (1996) Potent antisense oligonucleotides to the human multidrug resistance-1 mRNA are rationally selected by mapping RNA-accessible sites with oligonucleotide libraries. *Nucleic Acids Res*, 24: 1901–1907.
70. Scherr M and Rossi JJ. (1998) Rapid determination and quantitation of the accessibility to native RNAs by antisense oligodeoxynucleotides in murine cell extracts. *Nucleic Acids Res*, 26: 5079–5085.
71. Maher LJ and Dolnick BJ. (1988) Comparative hybrid arrest by tandem antisense oligodeoxyribonucleotides or oligodeoxyribonucleoside methylphosphonates in a cell-free system. *Nucleic Acids Res*, 16: 3341–3358.
72. Walter AE, Turner DH, et al. (1994) Coaxial stacking of helices enhances binding of oligoribonucleotides and improves predictions of RNA folding. *Proc Natl Acad Sci U S A*, 91: 9218–9222.
73. Buvoli A, Buvoli M, and Leinwand LA. (2000) Enhanced detection of tRNA isoacceptors by combinatorial oligonucleotide hybridization. *RNA*, 6: 912–918.
74. Krieg AM, Yi AK, and Hartmann G. (1999) Mechanisms and therapeutic applications of immune stimulatory CpG DNA. *Pharmacol Ther*, 84: 113–120.
75. Lipford GB, Bauer M, et al. (1997) CpG-containing synthetic oligonucleotides promote B and cytotoxic T cell responses to protein antigen: a new class of vaccine adjuvants. *Eur J Immunol*, 27: 2340–2344.
76. Sun S, Zhang X, et al. (1998) Type I interferon-mediated stimulation of T cells by CpG DNA. *J Exp Med*, 188: 2335–2342.
77. Bendigs S, Salzer U, et al. (1999) CpG-oligodeoxynucleotides co-stimulate primary T cells in the absence of antigen-presenting cells. *Eur J Immunol*, 29: 1209–1218.

78. Jakob T, Walker PS, et al. (1999) Bacterial DNA and CpG-containing oligodeoxynucleotides activate cutaneous dendritic cells and induce IL-12 production: implications for the augmentation of Th1 responses. *Int Arch Allergy Immunol*, 118: 457–461.
79. Krieg AM, Matson S, and Fisher E. (1996) Oligodeoxynucleotide modifications determine the magnitude of B cell stimulation by CpG motifs. *Antisense Nucleic Acid Drug Dev*. 6: 133–139.
80. Branda RF, Moore AL, et al. (1996) Amplification of antibody production by phosphorothioate oligodeoxynucleotides. *J Lab Clin Med*, 128: 329–338.
81. Benimetskaya L, Berton M, et al. (1997) Formation of a G-tetrad and higher order structures correlates with biological activity of the RelA (NF-kappaB p65) 'antisense' oligodeoxynucleotide. *Nucleic Acids Res*, 25: 2648–2656.
82. Wang W, Chen HJ, et al. (1998) A comparison of guanosine-quartet inhibitory effects versus cytidine homopolymer inhibitory effects on rat neointimal formation. *Antisense Nucleic Acid Drug Dev*, 8: 227–236.
83. Williamson JR, Raghuraman MK, and Cech TR. (1989) Monovalent cation-induced structure of telomeric DNA: the G-quartet model. *Cell*, 59: 871–880.
84. Schultze P, Macaya RF, and Feigon J. (1994) Three-dimensional solution structure of the thrombin-binding DNA aptamer d(GGTTGGTGTGGTTGG). *J Mol Biol*, 235: 1532–1547.
85. Chou SH, Zhu L, Reid BR. (1994) The unusual structure of the human centromere (GGA)<sub>2</sub> motif. Unpaired guanosine residues stacked between sheared G.A pairs. *J Mol Biol*, 244: 259–268.
86. Griffin LC, Toole JJ, and Leung LL. (1993) The discovery and characterization of a novel nucleotide-based thrombin inhibitor. *Gene*, 137: 25–31.
87. Wyatt JR, Vickers TA, et al. (1994) Combinatorially selected guanosine-quartet structure is a potent inhibitor of human immunodeficiency virus envelope-mediated cell fusion. *Proc Natl Acad Sci U S A*, 91: 1356–1360.
88. Tam RC, Lin CJ, et al. (1999) Inhibition of CD28 expression by oligonucleotide decoys to the regulatory element in exon 1 of the CD28 gene. *J Immunol*, 163: 4292–4299.
89. Bates PJ, Kahlon JB, et al. (1999) Antiproliferative activity of G-rich oligonucleotides correlates with protein binding. *J Biol Chem*, 274: 26369–26377.
90. Murchie AI and Lilley DM. (1994) Tetraplex folding of telomere sequences and the inclusion of adenine bases. *EMBO J*, 13: 993–1001.
91. Olivas WM and Maher LJ. (1995) Overcoming potassium-mediated triplex inhibition. *Nucleic Acids Res*, 23: 1936–1941.
92. Ratmeyer L, Vinayak R, et al. (1994) Sequence specific thermodynamic and structural properties for DNA.RNA duplexes. *Biochemistry*, 33: 5298–5304.
93. Tu GC, Cao QN, et al. (1998) Tetranucleotide GGGA motif in primary RNA transcripts. Novel target site for antisense design. *J Biol Chem*, 273: 25125–25131.